

MRI GROIN – SR ROUTINE

- Patient in supine position and place marker at site of pain
- Use large surface / body coil
- Coronal STIR of pelvis - include volume from anterior to the symphysis pubis through the entire sacrum
- Sagittals include volume centrally from medial acetabular walls
- Axials include volume from acetabular roofs distal to inferior pubic rami
- Obliques are parallel to anterior margin of iliac bone

SEQUENCE	FOV	Slice	TE	Matrix (min)
entire pelvis Cor STIR	38 cm	6 mm	40-60	256x256
Cor PD FS	16 cm	3.5 mm	20-30	256x256
Obl PD FS	16 cm	3.5 mm	20-30	256x256
Sag PD FS	16 cm	3.5 mm	20-30	256x256
Ax T1	16 cm	3.5 mm	min	256x256