

MRI PECTORALIS MAJOR– SR ROUTINE

► Use large FOV

Medial-lateral: cover to midline chest lateral: cover to midline chest *f*

Superior-inferior: cover down to mid-humeral humeral shaft

SEQUENCE	FOV	Slice	TE	Matrix (min)
Axial T1	30-36	3mm	min	256 x 256
Axial T2 FSE	30-36	3mm	80-100	256 x 256
Cor Obl STIR	30-36	3mm	20-40	256 x 256
Sag Oblique T2 FS FSE	30-36	3mm	40-60	256 x 256