

MRI SI JOINT/SACRUM– SR ROUTINE

- Patient in supine position as straight as possible
- Coronal obliques parallel to sacral body, cover bone only - include L3/L4 disc space at top of FOV to bottom of coccyx
- Axial obliques 90° to oblique coronal images - extend from L4/L5 disc space through end of coccyx
- No contrast unless otherwise specified or if for sacroiliitis

SEQUENCE	FOV	Slice	TE	Matrix (min)
Cor Obl T1	16 cm	3-4mm	min	320x256
Cor Obl STIR	16 cm	3-4mm	50	256x205
Ax Obl T1	16 cm	3-4mm	min	320x256
Ax Obl T2 FS	16 cm	3-4mm	80	256x205
Sag T1	16 cm	3-4mm	min	320x208
Sag STIR	16 cm	3-4 mm	50	256x208