

MRI STERNOCLAVICULAR JOINT– SR ROUTINE

- Patient in prone or supine position
- Surface coil
- Coil size depends on the examined area: sterno-clavicular joints

SEQUENCE	FOV	Slice	TE	Matrix (min)
Axial T1	250	3mm	min	400
Ax PD FS	250	3 mm	15-25	400
Cor PD FS	250	3mm	15-25	400
Sag T2 FS	250	3mm	40-60	400
CE T1 FS*	250	3	10-20	400

*if optimal assessment of synovitis/osteitis is needed